

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Isoprene - 99%, contains <1000 ppm p-tert-butylcatechol as inhibitor

Product Number: I19551
CAS Number: 78-79-5
MDL: MFCD00008600
Formula: C₅H₈
Formula Weight: 68.12 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Inhibitor	≤ 999 ppm
p-tert-Butylcatechol	

Specification: PRD.0.ZQ5.10000011390

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.