

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
NORIT(R) RX1.5 EXTRA

Product Number: 901939

Formula: C
Formula Weight: 12.01 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Powder or Crystals or Granules
Moisture Content	≤ 5 %(m)
Iron (Fe) (acid extr.)	≤ 300 mg/kg
Surface Area: BET	≥ 1370
Product of Supplier Cabot Corporation	Confirmed
Registered Trademark NORIT® is a registered trademark of Cabot Corporation	Note

Specification: PRD.0.ZQ5.10000094133

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.