

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
AW Standard Super-Cel® NFProduct Number: 221791  
CAS Number: 91053-39-3  
MDL: MFCD00132803

TEST	Specification
------	---------------

Appearance (Color)	Faint Beige to Beige and Faint Brown to Brown
Appearance (Form)	Powder
150 Mesh Tyler Sieve Retain	≤ 5.0 %
Loss on Drying	0.01 - 0.50 %
Registered Trademark	Confirmed
Super Cel is a registered trademark of Imerys Minerals California, Inc.	

Specification: PRD.0.ZQ5.10000043712

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.