

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Cyanopropan-2-yl N-methyl-N-(pyridin-4-yl)carbamodithioate - 97% (HPLC)

Product Number: **736236**
CAS Number: 1158958-96-3

Formula: C₁₁H₁₃N₃S₂
Formula Weight: 251.37 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Faint Yellow to Faint Brown
Appearance (Form)	Powder or Crystals
Proton NMR Spectrum	Conforms to Structure
Carbon	46.3 - 58.9 %
Nitrogen	14.7 - 18.7 %
Sulfur	22.4 - 28.6 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.2.ZQ5.10000013993

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.