

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Neodymium(III) chloride - anhydrous, powder, ≥99.99% trace metals basis

Product Number: 449946
CAS Number: 10024-93-8
MDL: MFCD00011132
Formula: Cl₃Nd
Formula Weight: 250.60 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Pale Green to Light Green	
Appearance (Form)	Powder
Complexometric EDTA	56.7 - 58.4 %
% Nd	
ICP Major Analysis	Confirmed
Confirms Neodymium Component	
Trace Rare Earth Analysis	≤ 100.0 ppm
Purity	Meets Requirements
≥99.99% Based On Rare Earth Analysis	

Specification: PRD.0.ZQ5.10000032640

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.