


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tyrphostin A9 - solid

Product Number: **T182**
CAS Number: 10537-47-0
MDL: MFCD00209853
Formula: C₁₈H₂₂N₂O
Formula Weight: 282.38 g/mol



TEST	Specification
Appearance (Color)	Light Yellow to Yellow
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Carbon	74.6 - 78.5 %
Nitrogen	9.7 - 10.2 %
Solubility (Color)	Light Yellow to Yellow
Solubility (Turbidity)	Clear
20 mg/mL, EtOH	
Purity (HPLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000054588

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.