


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methylphosphonic acid - 98%

Product Number: 289868
CAS Number: 993-13-5
MDL: MFCD00002137
Formula: CH₅O₃P
Formula Weight: 96.02 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Flake or Crystals or Crystalline Chunk(s)
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	97.5 - 102.5 %
GC (area %)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000024149

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.