


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tributyltin isocyanate

Product Number: 530786
CAS Number: 681-99-2
MDL: MFCD00019763
Formula: C₁₃H₂₇NOSn
Formula Weight: 332.07 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
ICP Major Analysis	Confirmed
Confirms Tin Component	

Specification: PRD.0.ZQ5.10000054382

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.