


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Pregnen-3 β -ol-20-one - $\geq 98\%$

Product Number: **P9129**
CAS Number: 145-13-1
MDL: MFCD00003628
Formula: C₂₁H₃₂O₂
Formula Weight: 316.48 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
at 10 mg/mL of ethanol	
Infrared spectrum	Conforms to Structure
Specific Rotation	27.0 - 30.0 °
c = 1 in ethanol, at 20 deg C	
Purity (HPLC)	$\geq 98\%$

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.