

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Alginate-RGD bioink – suitable for 3D bioprinting applications

Product Number: 901950

TEST	Specification
Syringe	-----
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Viscous Liquid
1H NMR Spectrum	Conforms to Structure
Bottle	-----
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Liquid
1H NMR Spectrum	Conforms to Structure
Expiration Date Period	-----
1 Year	

Specification: PRD.0.ZQ5.10000100210

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.