


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Phenolphthalein – pH indicator

Product Number: P9750
CAS Number: 77-09-8
MDL: MFCD00005913
Formula: C₂₀H₁₄O₄
Formula Weight: 318.32 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

50 mg/mL, EtOH

pH

Pass

8.0 (colorless) to 10 (red)

Specification: PRD.0.ZQ5.10000017831

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.