

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

2-Chloroethyl vinyl ether - 99%, contains triethanolamine as stabilizer, contains MEHQ as stabilizer

**Product Number:** 109983  
**CAS Number:** 110-75-8  
**MDL:** MFCD00000973  
**Formula:** C<sub>4</sub>H<sub>7</sub>ClO  
**Formula Weight:** 106.55 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Stabilizer	Confirmed
Contains 50-100ppm MEHQ + 300-500ppm Triethanolamine	

Specification: PRD.1.ZQ5.10000019888

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.