

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Potassium iodide – anhydrous, free-flowing, Redi-Dri™, ReagentPlus®, 99%

Product Number:

**793582**

CAS Number:

7681-11-0

Formula:

KI

Formula Weight:

166.01 g/mol

KI

### TEST

### Specification

Appearance (Color)

White

Appearance (Form)

Powder or Crystals

X-Ray Diffraction

Conforms to Structure

Titration by AgNO<sub>3</sub>

75.3 - 77.6 %

% I

Recommended Retest Period

42 Months

-----

Specification: PRD.1.ZQ5.10000049746

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.