

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Platinum(II) chloride - $\geq 99.9\%$

Product Number: 520632
CAS Number: 10025-65-7
MDL: MFCD00011181
Formula: Cl₂Pt
Formula Weight: 265.98 g/mol



TEST	Specification
Appearance (Color) Dark Green to Brown	Conforms to Requirements
Appearance (Form) Powder	Powder
Platinum (Pt) 71.5 - 73.5 %	71.5 - 73.5 %
ICP: Confirms Platinum Compone Confirmed	Confirmed
Trace Metal Analysis ≤ 1000.0 ppm	≤ 1000.0 ppm
Purity Meets Requirements $\geq 99.9\%$ Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.0.ZQ5.10000026245

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.