

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Buten-2-one - contains 0.3-1.0% hydroquinone as stabilizer, technical grade, 90%

Product Number:

M87509

CAS Number:

78-94-4

MDL:

MFCD00008777

Formula:


C₄H₆O

Formula Weight:

70.09 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Conforms to Requirements

Colorless to Dark Brown

Appearance (Form)

Liquid

Proton NMR spectrum

Conforms to Structure

Purity (GC)

≥ 88.0 %

Hydroquinone as Stabilizer

0.3 - 1.0 %

Water (by Karl Fischer)

≤ 7.0 %

Specification: PRD.2.ZQ5.10000004317

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.