

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

4-Hydroxy-2-(methylthio)-6-(trifluoromethyl)pyrimidine 97%

Product Number:

558818

CAS Number:

16097-62-4

MDL:


MFCD00153191

Formula:

C₆H₅F₃N₂OS

Formula Weight:

210.18 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Carbon

33.1- 35.5 %

Nitrogen

12.9- 13.8 %

Sulfur (S)

14.7- 15.8 %

Remarks:

Specification Date : 06.08.2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.