

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Piperidone monohydrate hydrochloride - 98%

Product Number:

151769

CAS Number:

40064-34-4

MDL:

MFCD00012776

Formula:


C₅H₉NO · HCl · H₂O

Formula Weight:

153.61 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

White to Yellow to Beige

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

97.5 - 102.5 %

Purity (TLC)

≥ 98 %

Solubility (Turbidity)

Clear

c = 25 mg/ml 1N HCl

Solubility (Color)

Colorless to Yellow

Specification: PRD.1.ZQ5.10000079525

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.