

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Poly(ethylene glycol) dimethacrylate – average  $M_n$  10,000, contains MEHQ as inhibitor

Product Number:

725684

CAS Number:

25852-47-5



Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

White to Yellow to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Solid

Proton NMR spectrum

Conforms to Structure

Loss on Drying

Molecular Number by GPC

8500 - 11500

(GPC, RI, THF, PEO)

Poly Dispersity MW/MN

 $\leq 1.1$ 

Specification: PRD.0.ZQ5.10000049540

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.