


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 437417
CAS Number: 35112-53-9
MDL: MFCD00040661
Formula: BaO3S2
Formula Weight: 249.46 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder to Crystalline Powder
ICP Major Analysis Confirms Barium & Sulfur Components	Confirmed
Titration by Iodine	97.5 - 102.5 %

Specification: PRD.0.ZQ5.10000048105

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.