

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Fibronectin solution from human fibroblasts – cell culture derived, ~0.5 mg/mL, sterile-filtered, BioReagent, suitable for cell culture

Product Number:**F0556****CAS Number:**

86088-83-7

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color) Colorless to faint yellow	Conforms
Appearance (Form)	Solution
Sterility	Sterile
Protein Content in CAPS saline buffer	0.4 - 0.6 mg/ml
SDS Electrophoresis Fibronectin band of 250 KDA detected	Conforms
Country of Origin	United States of America
Country of Manufacture	United States of America
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000012800

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.