


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Oil Red O solution - 0.5% in propylene glycol

Product Number: **O1516**



Formula Weight: 408.49 g/mol

| TEST | Specification |
|--------------------------|---------------|
| Appearance | Conforms |
| Dark Red, Sediment Free | |
| UV/VIS Absorbance | 0.40 - 1.00 |
| Lambda Max 506 to 520 nm | |
| Lambda Max | 506 - 520 nm |

Specification: PRD.0.ZQ5.10000058129

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.