

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium 1-(phenylsulfonyl)-1H-indole-3-trifluoroborate - 95%

Product Number:


750263

Formula:

C₁₄H₁₀BF₃KNO₂S

Formula Weight:

363.20 g/mol



TEST

Specification

Appearance (Color)

White to Yellow and Faint Beige to Beige
and Faint Brown to Brown

Appearance (Form)

Powder or Crystals

Proton NMR spectrum

Conforms to Structure

Fluorine NMR Spectra

Conforms to Structure

Boron NMR Spectra

Conforms to Structure

Carbon

43.5 - 49.1 %

Nitrogen

3.6 - 4.1 %

Sulfur

8.3 - 9.4 %

Specification: PRD.1.ZQ5.10000020759

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.