


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Thiabendazole - BioReagent, plant cell culture tested, powder

Product Number: T5535
CAS Number: 148-79-8
MDL: MFCD00005587
Formula: C₁₀H₇N₃S
Formula Weight: 201.25 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Faint Yellow to Brown
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/ml, Dimethylformide	
Purity (Titration by HCL04)	≥ 98.0 %
Purity (TLC)	≥ 99 %
Plant Cell Culture Test	Pass

Specification: PRD.0.ZQ5.10000047230

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.