

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bromine - reagent grade

Product Number: 207888
CAS Number: 7726-95-6
MDL: MFCD00010896
Formula: Br₂
Formula Weight: 159.81 g/mol



TEST	Specification
------	---------------

Appearance (Color)

Dark Red to Very Dark Red and Red-Brown
and Very Dark Red-Brown

Appearance (Form)

Liquid

Titration with Na₂S₂O₃

98.0 - 102.0 %

Specification: PRD.1.ZQ5.10000034514

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.