

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2xYT Medium - EZMix™ Powder microbial growth medium

Product Number: Y2627

TEST	Specification
------	---------------

Appearance (Color)	Light Yellow
Appearance (Form)	Powder
Granular powder	
Suitability	Confirmed
Greater than or equal to 80% of the growth observed in the control.	

Specification: PRD.2.ZQ5.1000009365

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.