

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethyl 1-bromocyclobutanecarboxylate - 96%

Product Number:**197297**

CAS Number:

35120-18-4

MDL:

MFCD00001322

Formula:


C₇H₁₁BrO₂

Formula Weight:

207.07 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Titration by AgNO₃ (BRCO)

95.5 - 104.5 %

% Purity Based On Bromide Content By Titration With AGNO₃ After O₂

Combustion

Purity (GC)

≥ 95.5 %

Specification: PRD.0.ZQ5.10000042627

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.