

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4,4'-(4,4'-Isopropylidenediphenoxy)bis(phthalic anhydride) - 97%

Product Number:

478032

CAS Number:

38103-06-9

MDL:


MFCD00319167

Formula:

C₃₁H₂₀O₈

Formula Weight:

520.49 g/mol



TEST

Specification

Appearance (Color)

Off White to Yellow

Appearance (Form)

Powder or Crystals or

Crystalline Chunk(s) or Granule or Flakes

Infrared Spectrum

Conforms to Structure

Carbon Content

70.0 - 75.1 %

Specification: PRD.0.ZQ5.10000100669

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.