


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Flutamide

Product Number: **F9397**
CAS Number: 13311-84-7
MDL: MFCD00072009
Formula: C₁₁H₁₁F₃N₂O₃
Formula Weight: 276.21 g/mol



TEST	Specification
Appearance (Color)	Yellow
Appearance (Form)	Powder
Solubility (Color)	Yellow to Yellow-Green
Solubility (Turbidity)	Clear to Hazy
50 mg/mL, EtOH	
Carbon	46.8 - 49.8 %
Nitrogen	9.8 - 10.4 %
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000043100

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.