

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Propylamine hydrochloride - ≥99%

Product Number: 242543
CAS Number: 556-53-6
MDL: MFCD00012911
Formula: C₃H₉N · HCl
Formula Weight: 95.57 g/mol

CH₃CH₂CH₂NH₂ · HCl

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Titration with AgNO ₃	99.0 - 101.0 %
Solubility (Turbidity) 100mg/ml, H ₂ O	Clear
Solubility (Color)	Colorless

Specification: PRD.2.ZQ5.10000004250

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.