

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
p-Xylene-bis(N-pyridinium bromide) - ≥95% (TLC)

Product Number: **D0940**
CAS Number: 14208-10-7

Formula: C₁₈H₁₈Br₂N₂
Formula Weight: 422.16 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
EmM	8.6 - 9.2
in H ₂ O	
Wavelength	256 - 262 nm
Lambda Max	
Purity (TLC)	≥ 95 %
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL, H ₂ O	

Specification: PRD.0.ZQ5.10000038914

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.