

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hydroxide – ultra dry, powder or crystals, 99.99% trace metals basis

Product Number: **757527**
CAS Number: 1310-73-2

NaOH

Formula Weight: 40.00 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
X-Ray Diffraction	Conforms to Structure
Titration by HCl	99.0 - 101.0 %
Purity	Meets Requirements
99.99% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.2.ZQ5.10000033011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.