


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Grubbs Catalyst® M204 - Umicore

Product Number: **569747**
 CAS Number: 246047-72-3
 MDL: MFCD03453237
 Formula: C₄₆H₆₅Cl₂N₂PRu
 Formula Weight: 848.97 g/mol
 Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color) Dark Red to Dark Brown	Conforms to Requirements
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	61.2 - 69.0 %
Nitrogen	3.1 - 3.5 %

Specification: PRD.1.ZQ5.10000009384

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

