

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum(0)-2,4,6,8-tetramethyl-2,4,6,8-tetravinylcyclotetrasiloxane complex solution - 0.104 M in methylvinylcyclsiloxanes

Product Number:

479543

CAS Number:

68585-32-0

MDL:


MFCD00166503

Formula:

C₁₂H₂₄O₄PtSi₄

Formula Weight:

539.74 g/mol



TEST

Specification

Appearance (Color)

Yellow to Brown

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

ICP Major Analysis

Confirmed

Confirms Platinum Component

Specification: PRD.0.ZQ5.1000020623

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.