

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Methyl-2-oxopentanoic acid sodium salt - ≥98%

Product Number: 198978
CAS Number: 3715-31-9
MDL: MFCD00002584
Formula: C₆H₉NaO₃
Formula Weight: 152.12 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder or Crystalline Powder	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	98.0 - 102.0 %
Solubility (Turbidity) 50mg/ml, H ₂ O	Clear
Solubility (Color)	Colorless
Purity (TLC)	≥ 98 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.