

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: **F1168**

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from a 0.2 um filtered solution of PBS containing 50 ug BSA/ug FGF-9	
Bioactivity	0.1 - 5.0 ng/ml
ED50 (cell proliferation, BALB/3T3 cells)	
Endotoxin Assay	Conforms
< 1.0 EU/ug FGF-9	
Purity by SDS Electrophoresis	Conforms
> 97%	

Specification: PRD.0.ZQ5.10000101804

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.