


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Pyridinecarboxylic acid - 97%

Product Number: 297763
CAS Number: 50681-25-9
MDL: MFCD00011576
Formula: C₅H₄N₂O₂
Formula Weight: 124.10 g/mol



TEST

Specification

Appearance (Color)	White to Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 96.5 %
Titration with NaOH	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000070625

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.