

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hyaluronidase from sheep testes - Type V, lyophilized powder, $\geq 1,500$ units/mg solid

Product Number:

H6254

CAS Number:

37326-33-3

MDL:

MFCD00131351

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/mg solid

 ≥ 1500

Unit Definition: One unit will cause a change in % transmittance at 600 nm of 0.330 per minute at pH 5.35 at 37 deg C in a 2.0 mL reaction mixture (45 minute assay).

Specification: PRD.2.ZQ5.10000026633

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.