

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Collagen, Type I solution from rat tail - BioReagent, suitable for cell culture, sterile-filtered

 Product Number:
 C3867

 CAS Number:
 9007-34-5

 MDL:
 MFCD00130825

WII 0500 10002

Storage Temperature: 2 - 8 ℃

TEST Specification

Appearance (Color) Colorless

Appearance (Turbidity)

Clear to Slightly Hazy

Appearance (Form) Liquid

at ambient temperature

Sterility (Direct Inoculation) Conforms

by USP.

Negative for the presence of bacteria and fungi.

Mycoplasma Test by Barile None Detected Protein by Biuret  $$\geq 0.5~{\rm mg/ml}$$ 

SDS Electrophoresis Consistent with Past Lots

SDS-Page

Expiration Date Period ------

1 year

Specification: PRD.2.ZQ5.10000007975

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.