

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hydrobromic acid - 48 wt. % in H₂O, ≥99.99%

Product Number:

339245

CAS Number:

10035-10-6

MDL:

MFCD00011323

Formula:

HBr

Formula Weight:

80.91 g/mol

HBr**TEST****Specification**

Appearance (Color)

Colorless to Very Dark Yellow to Very Dark Brown-Yellow

Appearance (Form)

Liquid

Purity (Titration by NaOH)

47.0 - 49.0 %

Purity

Meets Requirements

> = 99.99% Based On Trace Metals Analysis

Trace Metal Analysis

≤ 100.0 ppm

Specification: PRD.1.ZQ5.10000028565

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.