

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
4-Nitrophenol solution - 10 mM

Product Number: **N7660**

MDL: MFCD00007331  
Formula: C<sub>6</sub>H<sub>5</sub>NO<sub>3</sub>  
Formula Weight: 139.11 g/mol  
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Essentially free of particles	
pH:	3.5 - 3.7
1X at 25°C	
Absorbance	0.890 - 0.924
Mean Test A401	

Specification: PRD.1.ZQ5.10000043869

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.