

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Gold etchant, nickel compatible

Product Number: 651842

MDL: MFCD07370792

TEST Specification

Appearance (Color)
Appearance (Form)
ICP Major Analysis
Confirms P and K Components
Titration by HCL
Wt.% (NH4)2HPO4
Titration by Na2S2O3
Wt.% I2
Titration by AgNO3
Wt.% KI

Specification: PRD.2.ZQ5.10000038277

Yellow to Dark Yellow to Dark Orange

Liquid Confirmed

4 - 10 %

5 - 10 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.