

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Copper(II) sulfate - anhydrous, powder, ≥99.99% trace metals basis

Product Number: 451657 CAS Number: 7758-98-7

MDL: MFCD00010981

Formula: CuO4S

Formula Weight: 159.61 g/mol

TEST Specification

Appearance (Color) Conforms to Requirements

White to Light Grey Appearance (Form)

ICP Major Analysis

Confirms Copper Component

Trace Metal Analysis

Purity

> = 99.99% Based On Trace Metals Analysis

Specification: PRD.0.ZQ5.10000023820

< 100.0 ppm

Pow der

Confirmed

CuSO₄

Meets Requirements

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.