

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Phosphate buffered saline – tablet

Product Number: P4417

MDL: MFCD00131855

TEST	Specification
Total Tablet Weight	1814.5 - 2005.5 mg/tab
Appearance	White
Form	Round Tablet
Diameter	5/8 inches
Dissolution Time	≤ 10.00 min
pH	7.20 - 7.60
Note	-----

Dissolution Time and pH test results from a solution of one tablet in 200 ml of purified water.

Specification: PRD.2.ZQ5.10000006798.000

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.