

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Carmine - certified by the Biological Stain Commission

Product Number: **C6152**
CAS Number: 1390-65-4

TEST	Specification
Appearance (Color)	Red
Appearance (Form)	Powder
Solubility (Color) 1 mg/mL in 2 N NH ₄ OH	Red
Infrared spectrum	Conforms to Structure
Biological Stain Commission	Pass

Specification: PRD.0.ZQ5.10000028521

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.