

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium phosphate dibasic – BioReagent, suitable for cell culture, suitable for insect cell culture, ≥99.0%

Product Number:

S5136

CAS Number:

7558-79-4

MDL:


MFCD00003496

Formula:

HNa₂O₄P

Formula Weight:

141.96 g/mol



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H ₂ O	
Loss on Drying	≤ 1.0 %
Titration with NaOH	≥ 99.0 %
Dry Basis	
Cell Culture Test	Pass
Insect Cell Test	Pass
Recommended Retest Period	-----
3 years	

Specification: PRD.1.ZQ5.10000011402

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.