

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trypsin from porcine pancreas - lyophilized powder, BioReagent, suitable for cell culture, 1,000-2,000 BAE units/mg solid

Product Number: T4799
CAS Number: 9002-07-7
MDL: MFCD00082094

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow to Light Brown
Appearance (Form)	Powder
Solubility (Turbidity)	Clear
35 mg/ml, HBSS	
BAEE units/mg solid	1000 - 2000
BTEE units/mg solid	≤ 10
Cell Dissociation Assay	Pass
Mycoplasma Test by Barile	Negative
Porcine Parvovirus	None Detected
by 9CFR	
Recommended Retest Period	-----
4 Years	

Specification: PRD.1.ZQ5.10000035684

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.