


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Acryloyl chloride - 97.0%, contains <210 ppm MEHQ as stabilizer

Product Number: 549797
CAS Number: 814-68-6
MDL: MFCD00000731
Formula: C₃H₃ClO
Formula Weight: 90.51 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Pale Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 97.0 %
Stabilizer	≤ 210 ppm
MEHQ	
Recommended Retest Period	-----
1 Year	

Specification: PRD.2.ZQ5.10000008148

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.