

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sulfur trioxide – contains stabilizer, ≥99%

Product Number:**425478**

CAS Number:

7446-11-9

MDL:


MFCD00011452

Formula:

O₃S

Formula Weight:

80.06 g/mol



TEST	Specification
------	---------------

Supplier - Purity

≥ 99.0 %

Sulfur Trioxide

Stabilizer

Confirmed

Note

Note

Appearance may be slightly turbid
with a straw or tan colored tint

Specification: PRD.8.ZQ5.10000011889

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.