


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cesium acetate - Vetec[™] reagent grade

Product Number: V900146
CAS Number: 3396-11-0
Formula: C₂H₃CsO₂
Formula Weight: 191.95 g/mol



TEST

Specification

Appearance (Colour)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Solubility (Colour)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Very Slightly Hazy
C = 20%, H ₂ O	
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	95.0 - 105.0 %

Specification: PRE.0.ZQ5.10000002617

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.