

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Indium tin oxide coated PET – surface resistivity 200 Ω/sq, L × W × thickness 1 ft × 1 ft × 5 mil, sheet

Product Number:

749745



TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form) Sheet	Confirmed
Size 1ft x 1ft x 5 mil	Confirmed
Coating In ₂ O ₃ /(SnO ₂) _x	Confirmed
Surface Resistivity (Ohm.sq) Specification: Report Result	
Transmittance at 550nm Specification: Report Result	

Specification: PRD.0.ZQ5.10000018271

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.